

## Product Datasheet

### Human PDE3B ELISA Kit (orb551897)

|                          |                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Phosphodiesterase 3B, cGMP<br>Inhibited (PDE3B) Sandwich ELISA Kit                                      |
| <b>Reactivity</b>        | Human                                                                                                         |
| <b>Range</b>             | 78.125-5000pg/mL                                                                                              |
| <b>Form/Appearance</b>   | Traditional                                                                                                   |
| <b>Storage</b>           | All reagents should be kept according to the<br>instructions on the packaging upon receipt.                   |
| <b>Note</b>              | For research use only                                                                                         |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell<br>lysates, cell culture supernates or other<br>biological fluids.    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell<br>lysates, cell culture supernatants and other<br>biological fluids. |
| <b>Sensitivity</b>       | 26.35pg/mL                                                                                                    |
| <b>Expiration Date</b>   | Please enquire.                                                                                               |