

Product Datasheet

Human PFKM ELISA Kit (Ready to Use) (orb551901)

Description	Human Phosphofructokinase, Muscle (PFKM) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P08237
Sensitivity	0.114ng/mL
Entrez	5213
Expiration Date	Please enquire.