

## Product Datasheet

### Human PDGFAA ELISA Kit (Ready to Use) (orb551952)

|                          |                                                                                                                                                                       |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Platelet Derived Growth Factor<br>AA (PDGFAA) Sandwich ELISA KitReady<br>to Use: Pre-diluted detecti...                                                         |
| <b>Reactivity</b>        | Human                                                                                                                                                                 |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                                                                         |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                                                                          |
| <b>Storage</b>           | All reagents should be stored at -20°C<br>upon receipt                                                                                                                |
| <b>Note</b>              | For research use only                                                                                                                                                 |
| <b>Application notes</b> | Pre-diluted detection reagents and<br>enhanced stabilityserum, plasma,<br>tissue homogenates, cell lysates, cell<br>culture supernates or other biological<br>fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates,<br>cell lysates, cell culture supernatants<br>and other biological fluids.                                                         |
| <b>Sensitivity</b>       | 0.055ng/mL                                                                                                                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                                                                                                                       |