

## Product Datasheet

### Human PAWP ELISA Kit (Ready to Use) (orb551986)

|                          |                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Postacrosomal Sheath WW<br>Domain Binding Protein (PAWP)<br>Sandwich ELISA KitReady to Use: Pre-<br>dilu... |
| <b>Reactivity</b>        | Human                                                                                                             |
| <b>Range</b>             | 31.2-2000pg/mL                                                                                                    |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                      |
| <b>Storage</b>           | All reagents should be stored at -20°C<br>upon receipt                                                            |
| <b>Note</b>              | For research use only                                                                                             |
| <b>Application notes</b> | Pre-diluted detection reagents and<br>enhanced stabilitytissue homogenates<br>or other biological fluids.         |
| <b>Sample Types</b>      | tissue homogenates and other<br>biological fluids.                                                                |
| <b>Sensitivity</b>       | 12.4pg/mL                                                                                                         |
| <b>Expiration Date</b>   | Please enquire.                                                                                                   |