

Product Datasheet

Human PICP ELISA Kit (Ready to Use) (orb552012)

Description	Human Procollagen I C-Terminal Propeptide (PICP) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	23.1pg/mL
Expiration Date	Please enquire.