

Product Datasheet

Human PLOD3 ELISA Kit (Ready to Use) (orb552023)

Description	Human Procollagen Lysine-2-Oxoglutarate-5-Dioxygenase 3 (PLOD3) Sandwich ELISA KitReady to Use: Pre-...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O60568
Sensitivity	0.106ng/mL
Entrez	8985
Expiration Date	Please enquire.