

Product Datasheet

Human PTPLA ELISA Kit (Ready to Use) (orb552122)

Description	Human Protein Tyrosine Phosphatase Like Protein A (PTPLA) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	B0YJ81
Sensitivity	0.25ng/mL
Entrez	9200
Expiration Date	Please enquire.