

Product Datasheet

Human PTP4A3 ELISA Kit (Ready to Use) (orb552126)

Description	Human Protein Tyrosine Phosphatase Type IVA 3 (PTP4A3) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	075365
Sensitivity	0.118ng/mL
Entrez	11156
Expiration Date	Please enquire.