

Product Datasheet

Human RAB1A ELISA Kit (Ready to Use) (orb552160)

Description	Human RAB1A, Member RAS Oncogene Family (RAB1A) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P62820
Sensitivity	0.067ng/mL
Entrez	5861
Expiration Date	Please enquire.