

Product Datasheet

Human RORg ELISA Kit (Ready to Use) (orb552175)

Description	Human RAR Related Orphan Receptor Gamma (RORg) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P51449
Sensitivity	0.071ng/mL
Entrez	6097
Expiration Date	Please enquire.