

Product Datasheet

Human RNASEH ELISA Kit (Ready to Use) (orb552240)

Description	Human Ribonuclease H (RNASEH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O60930
Sensitivity	0.052ng/mL
Entrez	246243
Expiration Date	Please enquire.