

## Product Datasheet

### Human SSH2 ELISA Kit (Ready to Use) (orb552398)

|                          |                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Slingshot Homolog 2 (SSH2)<br>Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ... |
| <b>Reactivity</b>        | Human                                                                                                      |
| <b>Range</b>             | 0.312-20ng/mL                                                                                              |
| <b>Form/Appearance</b>   | Ready to Use                                                                                               |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                        |
| <b>Note</b>              | For research use only                                                                                      |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.        |
| <b>Sample Types</b>      | tissue homogenates and other biological fluids.                                                            |
| <b>Uniprot ID</b>        | <b>Q76I76</b>                                                                                              |
| <b>Sensitivity</b>       | 0.113ng/mL                                                                                                 |
| <b>Entrez</b>            | <b>85464</b>                                                                                               |
| <b>Expiration Date</b>   | Please enquire.                                                                                            |