

Product Datasheet

Human SGLT1 ELISA Kit (Ready to Use) (orb552411)

Description	Human Sodium/Glucose Cotransporter 1 (SGLT1) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P13866
Sensitivity	0.056ng/mL
Entrez	6523
Expiration Date	Please enquire.