

Product Datasheet

PI 3 Kinase Class 3 antibody [Out of stock] (orb55245)

Description

Rabbit polyclonal antibody against PI 3

Species/Host	Rabbit
Conjugation	PerCP
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide from the N-terminal of Human PI 3 Kinase Class 3
Target	PI 3 Kinase Class 3
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-661 was sufficient for detection of PI 3 Kinase Class 3 in 15 µg of mouse liver cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	101.5 kDa
Uniprot ID	Q8NEB9
NCBI	NP_002638.2
Entrez	5289
Dilution Range	WB (1:1000), ICC/IF (1:100); IHC (1:50)
Expiration Date	12 months from date of receipt.