

Product Datasheet

Human SGPP2 ELISA Kit (orb552450)

Description	Human Sphingosine-1-Phosphate Phosphatase 2 (SGPP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.059ng/mL
Expiration Date	Please enquire.