

Product Datasheet

Human TERF2IP ELISA Kit (Ready to Use) (orb552549)

Description	Human Telomeric Repeat Binding Factor 2 Interacting Protein (TERF2IP) Sandwich ELISA KitReady to Use...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9NYB0
Sensitivity	0.66ng/mL
Entrez	54386
Expiration Date	Please enquire.