

Product Datasheet

Human TLL1 ELISA Kit (Ready to Use) (orb552623)

Description	Human Tolloid Like Protein 1 (TLL1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O43897
Sensitivity	0.216ng/mL
Entrez	7092
Expiration Date	Please enquire.