

Product Datasheet

Human TRPM1 ELISA Kit (Ready to Use) (orb552651)

Description	Human Transient Receptor Potential Cation Channel Subfamily M, Member 1 (TRPM1) Sandwich ELISA KitRe...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q7Z4N2
Sensitivity	0.051ng/mL
Entrez	4308
Expiration Date	Please enquire.