

Product Datasheet

Human TPP2 ELISA Kit (Ready to Use) (orb552679)

Description	Human Tripeptidyl-Peptidase 2 (TPP2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.130ng/mL
Expiration Date	Please enquire.