

Product Datasheet

Human TP53BP1 ELISA Kit (Ready to Use) (orb552708)

Description	Human Tumor Protein p53 Binding Protein 1 (TP53BP1) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q12888
Sensitivity	0.138ng/mL
Entrez	7158
Expiration Date	Please enquire.