

Product Datasheet

Human UBE2L3 ELISA Kit (Ready to Use) (orb552731)

Description	Human Ubiquitin Conjugating Enzyme E2L3 (UBE2L3) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P68036
Sensitivity	0.072ng/mL
Entrez	7332
Expiration Date	Please enquire.