

Product Datasheet

Human UGT2B7 ELISA Kit (Ready to Use) (orb552740)

Description	Human UDP Glucuronosyltransferase 2 Family, Polypeptide B7 (UGT2B7) Sandwich ELISA KitReady to Use: ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P16662
Sensitivity	0.113ng/mL
Entrez	7364
Expiration Date	Please enquire.