

Product Datasheet

Human UGGT1 ELISA Kit (Ready to Use) (orb552742)

Description	Human UDP-Glucose Glycoprotein Glucosyltransferase 1 (UGGT1) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NYU2
Sensitivity	0.113ng/mL
Entrez	56886
Expiration Date	Please enquire.