

Product Datasheet

Human VIPR1 ELISA Kit (Ready to Use) (orb552787)

Description	Human Vasoactive Intestinal Peptide Receptor 1 (VIPR1) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P32241
Sensitivity	0.045ng/mL
Entrez	7433
Expiration Date	Please enquire.