

Product Datasheet

Human VACHT ELISA Kit (Ready to Use) (orb552797)

Description	Human Vesicular Acetylcholine Transporter (VACHT) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q16572
Sensitivity	0.117ng/mL
Entrez	6572
Expiration Date	Please enquire.