


Product Datasheet

IKK gamma (phospho-Ser31) antibody (orb5528)

Description	IKK gamma (phospho-Ser31) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IKK gamma around the phosphorylation site of Ser31 (EE(p-S)PL)
Target	phospho-IKK gamma (Ser31)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	Q9Y6K9
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
pancreatic
carci...