

Product Datasheet

Human WFDC5 ELISA Kit (Ready to Use) (orb552819)

Description	Human WAP Four Disulfide Core Domain Protein 5 (WFDC5) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8TCV5
Sensitivity	35pg/mL
Entrez	149708
Expiration Date	Please enquire.