

Product Datasheet

Human XRCC6 ELISA Kit (Ready to Use) (orb552845)

Description	Human X-Ray Repair Cross Complementing 6 (XRCC6) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P12956
Sensitivity	0.112ng/mL
Entrez	2547
Expiration Date	Please enquire.