

Product Datasheet

Monkey IL1R2 ELISA Kit (Ready to Use) (orb552863)

Description	Monkey Interleukin 1 Receptor Type II (IL1R2) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Monkey
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8HXR9
Sensitivity	0.213ng/mL
Entrez	711638
Expiration Date	Please enquire.