

Product Datasheet

Mouse ADK ELISA Kit (Ready to Use) (orb552891)

Description	Mouse Adenosine Kinase (ADK) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P55264
Sensitivity	0.054ng/mL
Entrez	11534
Expiration Date	Please enquire.