

Product Datasheet

Mouse ADCY10 ELISA Kit (Ready to Use) (orb552892)

Description	Mouse Adenylate Cyclase 10, Soluble (ADCY10) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8C0T9
Sensitivity	0.132ng/mL
Entrez	271639
Expiration Date	Please enquire.