

Product Datasheet

Mouse beta 1 Adrenergic Receptor ELISA Kit (Ready to Use) (orb552896)

Description	Mouse Adrenergic Receptor Beta 1 (beta 1 Adrenergic Receptor) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P34971
Sensitivity	0.059ng/mL
Entrez	11554
Expiration Date	Please enquire.