

Product Datasheet

Mouse ALDH1A1 ELISA Kit (Ready to Use) (orb552908)

Description	Mouse Aldehyde Dehydrogenase 1 Family, Member A1 (ALDH1A1) Sandwich ELISA KitReady to Use: Pre- dilut...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P24549
Sensitivity	0.059ng/mL
Entrez	11668
Expiration Date	Please enquire.