

Product Datasheet

Mouse ALDM ELISA Kit (Ready to Use) (orb552911)

Description	Mouse Aldehyde Dehydrogenase, Mitochondrial (ALDM) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P47738
Sensitivity	0.12ng/mL
Entrez	11669
Expiration Date	Please enquire.