

Product Datasheet

Mouse ALDOC ELISA Kit (Ready to Use) (orb552915)

Description	Mouse Aldolase C, Fructose Bisphosphate (ALDOC) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.278ng/mL
Expiration Date	Please enquire.