

Product Datasheet

Mouse a1AGP ELISA Kit (Ready to Use) (orb552921)

Description	Mouse Alpha-1-Acid Glycoprotein (a1AGP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Mouse
Range	62.5-4000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	J3JRU4
Sensitivity	19.3ng/mL
Entrez	18407
Expiration Date	Please enquire.