

Product Datasheet

Mouse aZGP1 ELISA Kit (Ready to Use) (orb552924)

Description	Mouse Alpha-2-Glycoprotein 1, Zinc Binding (aZGP1) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	1.25-80ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q64726
Sensitivity	0.48ng/mL
Entrez	12007
Expiration Date	Please enquire.