

Product Datasheet

Mouse AIMP1 ELISA Kit (Ready to Use) (orb552934)

Description	Mouse Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) Sandwich ELISA...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q3UZG4
Sensitivity	14.8pg/mL
Entrez	13722
Expiration Date	Please enquire.