

Product Datasheet

Mouse APBB3 ELISA Kit (Ready to Use) (orb552944)

Description	Mouse Amyloid Beta Precursor Protein Binding Protein B3 (APBB3) Sandwich ELISA KitReady to Use: Pre-...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.255ng/mL
Expiration Date	Please enquire.