

Product Datasheet

Mouse AGGF1 ELISA Kit (Ready to Use) (orb552947)

Description	Mouse Angiogenic Factor With G Patch And FHA Domains 1 (AGGF1) Sandwich ELISA KitReady to Use: Pre-d...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q7TN31
Sensitivity	34pg/mL
Entrez	66549
Expiration Date	Please enquire.