

Product Datasheet

Mouse Angiopoietin 2 ELISA Kit (Ready to Use) (orb552950)

Description	Mouse Angiopoietin 2 (Angiopoietin 2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	035608
Sensitivity	6.1pg/mL
Entrez	11601
Expiration Date	Please enquire.