

Product Datasheet

Mouse ASNS ELISA Kit (Ready to Use) (orb553002)

Description	Mouse Asparagine Synthetase (ASNS) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q61024
Sensitivity	0.108ng/mL
Entrez	27053
Expiration Date	Please enquire.