

Product Datasheet

Mouse AST2 ELISA Kit (Ready to Use) (orb553004)

Description	Mouse Aspartate Aminotransferase 2 (AST2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P05202
Sensitivity	0.29ng/mL
Entrez	14719
Expiration Date	Please enquire.