

Product Datasheet

Mouse ABCC6 ELISA Kit (Ready to Use) (orb553008)

Description	Mouse ATP Binding Cassette Transporter C6 (ABCC6) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9R1S7
Sensitivity	0.124ng/mL
Entrez	27421
Expiration Date	Please enquire.