

## Product Datasheet

### Mouse ABCG2 ELISA Kit (Ready to Use) (orb553009)

|                          |                                                                                                                  |
|--------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse ATP Binding Cassette<br>Transporter G2 (ABCG2) Sandwich<br>ELISA KitReady to Use: Pre-diluted<br>detect... |
| <b>Reactivity</b>        | Mouse                                                                                                            |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                    |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                     |
| <b>Storage</b>           | All reagents should be stored at -20°C<br>upon receipt                                                           |
| <b>Note</b>              | For research use only                                                                                            |
| <b>Application notes</b> | Pre-diluted detection reagents and<br>enhanced stabilitytissue homogenates<br>or other biological fluids.        |
| <b>Sample Types</b>      | tissue homogenates and other<br>biological fluids.                                                               |
| <b>Uniprot ID</b>        | <b>Q7TMS5</b>                                                                                                    |
| <b>Sensitivity</b>       | 0.055ng/mL                                                                                                       |
| <b>Entrez</b>            | <b>26357</b>                                                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                                                  |