

Product Datasheet

Mouse ATP2A2 ELISA Kit (Ready to Use) (orb553010)

Description	Mouse ATPase, Ca++ Transporting, Cardiac Muscle, Slow Twitch 2 (ATP2A2) Sandwich ELISA KitReady to U...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O55143
Sensitivity	0.059ng/mL
Entrez	11938
Expiration Date	Please enquire.