

Product Datasheet

Mouse ATP6AP2 ELISA Kit (orb553013)

Description	Mouse ATPase, H ⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.055ng/mL
Expiration Date	Please enquire.