

Product Datasheet

Mouse ATP6AP2 ELISA Kit (Ready to Use) (orb553014)

Description	Mouse ATPase, H ⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) Sandwich ELISA KitReady to Us...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.055ng/mL
Expiration Date	Please enquire.