

Product Datasheet

Mouse ATP1a2 ELISA Kit (Ready to Use) (orb553019)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Alpha 2 Polypeptide (ATP1a2) Sandwich ELISA KitReady to Use: Pre-d...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q6PIE5
Sensitivity	0.042ng/mL
Entrez	98660
Expiration Date	Please enquire.