

Product Datasheet

Mouse ATP1b3 ELISA Kit (Ready to Use) (orb553021)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Beta 3 Polypeptide (ATP1b3) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P97370
Sensitivity	0.066ng/mL
Entrez	11933
Expiration Date	Please enquire.